



AN **EASY**, EFFICIENT AND AFFORDABLE SHOPRIDER MOBILITY DEVICE MEANS **FREEDOM** by Shoprider



Model Avanti

This **new compact model** combining the features of the power wheelchair and the scooter is suitable for travelling in urban areas or indoors such as art galleries, museums, airports, train stations, shopping malls etc. The standard version is equipped with SLA battery, yet an option with environment friendly battery technology - Lithium iron phosphate battery LiFePO4 - is also available.

SCOOCHAIR™ MINI (Model SC1331)	Technical specifications	
Overall Dimensions: L ¹ x W x H	mm	850 x 500 x 850
Number, Size of Tires	Front	2 x 6" Punch Proof Tires
	Rear	2 x 7" Punch Proof Tires
Total Weight (excluding Battery)		23,5
Total Weight (including Battery)	kg	31
Max. User Weight (on the level road)	kg	90
Battery Capacity	_V_Ah x Pcs	SLA 12V 10Ah x 2
Charger	Amp	2A (Off-Board)
Max. Speed* ² (on the level road)	km/h	Approx. 5~6
Safe Climbing Angle	°(%)	6° (10%)
Range* ³	km	14~16
Colors available	-	Red, Blue
Turning Radius	mm	835




31 kg



16 km/h

SHOPRIDER reserves all rights to change the wording and specification without prior notice. Data shown above is based on the time when tested.

*1. Include the rear pull wheel. *2. Driver weight may exceed weight of the unit, so speed must be reduced when turning. *3. Range was obtained and tested with 70kg weight load including tester & dummy. The actual driving range may vary by the factors shown below. (1) the weight of occupant; (2) level ground/gradient; (3) battery conditions; (4) type of charger; (5) ambient temperature; (6) the way of driving; (7) terrain; (8) etc.